

# Formalin Safety Just Got Safer



**Tissue-Tek® FormaGO®**

**Formalin Safety System**

# Tissue-Tek® FormaGO®

# The Formalin Safety System

The Tissue-Tek® FormaGO® Formalin Safety System includes products that further protect your staff and the environment from the hazardous effects of formalin.

Sakura developed this new product line to improve performance, reliability and convenience in handling formalin, with the safety of our valued customers and the environment in mind.

## USER SAFETY

### Tissue-Tek® FormaGO® Formalin Absorption Pads

- Sturdy and durable grossing pads with liner protect laboratory personnel by absorbing 18% more formalin during grossing<sup>1</sup>
- Now available in 3 sizes (small, medium and large) to help prevent specimen cross-contamination between cases



### Tissue-Tek® FormaGO® Formalin Absorption Wipes

- Ideal for small formalin spills
- Protects laboratory personnel from exposure to dangerous vapors
- Thicker, sturdier cloth-like wipes for better handling<sup>2</sup>



### Tissue-Tek® FormaGO® Formalin Absorption Granules

- Ideal for large formalin spills
- Protects laboratory personnel from exposure to dangerous vapors
- 35% lighter in weight<sup>3</sup> for more convenient handling
- Based on an innovative formulation, the granules retain their form making it easier to pick them up with a broom or brush, while leaving less residue after cleaning the spill<sup>3</sup>



## ENVIRONMENTAL SAFETY

### Tissue-Tek® FormaGO® Formalin Neutralizer

- Dissolves and neutralizes formalin up to 2 to 3 times faster - as low as in 5 minutes<sup>4</sup>
- Each premeasured pouch neutralizes 1 gallon of 10% formalin
- Dissolves with less physical effort<sup>4</sup>



### Tissue-Tek® FormaGO® Formaldehyde Analysis Kit

- Eliminates upfront pH measurement<sup>5</sup>
- Long 2-year shelf life improves inventory management
- Consistent results



### Tissue-Tek® FormaGO® Formalin Neutralizer Container

- Narrow opening reduces exposure to fumes
- Comfortable handle makes it convenient to carry



## COSTS

Now, more cost-effective than ever before

# Tissue-Tek® FormaGO® The Formalin Safety System

The Tissue-Tek® FormaGO® Formalin Safety System protects laboratory personnel and the environment by safely absorbing/neutralizing formalin.

The new Tissue-Tek® FormaGO® Formalin Safety System is safe, convenient and cost effective:

- 18% more absorption of formalin<sup>1</sup>
- 2 to 3 times faster neutralization of formalin - as low as in 5 minutes<sup>4</sup>
- Multiple sizes of pads to help prevent specimen cross-contamination between cases

*Time to switch to Tissue-Tek® FormaGO®!*

## NEW Tissue-Tek® FormaGO® products

| New Products   | New Product Codes |
|--|-------------------|
| Tissue-Tek® FormaGO® Formalin Starter Kit  | 9150              |
| Tissue-Tek® FormaGO® Formalin Neutralizer, (16 pouches)                              | 9151              |
| Tissue-Tek® FormaGO® Formaldehyde Analysis Kit<br>(2 reagent bottles and 100 strips) | 9152              |
| Tissue-Tek® FormaGO® Formalin Absorption Wipes,<br>(4 packs of 25 wipes)             | 9154              |
| Tissue-Tek® FormaGO® Formalin Absorption Pads, Small<br>(10 bags of 100 pads)        | 9155              |
| Tissue-Tek® FormaGO® Formalin Absorption Pads, Medium<br>(10 bags of 20 pads)        | 9156              |
| Tissue-Tek® FormaGO® Formalin Absorption Pads, Large<br>(10 bags of 10 pads)         | 9157              |
| Tissue-Tek® FormaGO® Formalin Absorption Granules<br>(2 jars)                        | 9158              |
| Tissue-Tek® FormaGO® Formalin Neutralizer Container (9.5L)                           | 9159              |

## Tissue-Tek® FormaGO® products replace these Tissue-Tek® Neutra-Guard® products

| Products                                     | (Product Codes) |
|--|-----------------|
| Tissue-Tek® Neutralex®<br>Starter Pack       | (4046)          |
| Tissue-Tek® Neutralex®                       | (4047)          |
| Tissue-Tek® Neutralex®<br>Aldehyde Test Kit  | (4050)          |
| Tissue-Tek® Neutra-Wipes®                    | (4279)          |
| Tissue-Tek® Neutra-Pads®                     | (4295)          |
| Tissue-Tek® Neutra-Form®                     | (4304)          |
| Tissue-Tek® Neutralex®<br>Reaction Container | (4276)          |

1 Comparison of Tissue-Tek® FormaGO® Formalin Absorption Pads Large with Tissue-Tek® Neutra-Pads®. 2 Comparison of Tissue-Tek® FormaGO® Formalin Absorption Wipes with Tissue-Tek® Neutra-Wipes®. 3 Comparison of Tissue-Tek® FormaGO® Formalin Absorption Granules with Tissue-Tek® Neutra-Form®. 4 Comparison of Tissue-Tek® FormaGO® Formalin Neutralizer with Tissue-Tek® Neutralex®. 5 Comparison of Tissue-Tek® FormaGO® Formaldehyde Analysis Kit with Tissue-Tek® Neutralex® Aldehyde Test Kit. Neutralex, Neutra-Wipes, Neutra-Pads, Neutra-Form and Neutra-Guard are trademarks of Scigen Inc.. Sakura, Tissue-Tek and FormaGO are trademarks of Sakura Finetek U.S.A., Inc.



CONTINUOUS INNOVATION FOR PATHOLOGY



Call to speak with a representative today at **(800) 725-8723**  
or visit [www.sakura-america.com](http://www.sakura-america.com)

Sakura Finetek USA, Inc.  
1750 West 214<sup>th</sup> Street  
Torrance, CA 90501 U.S.A.

0007792-01 RevA